


Product Name: Recombinant Human TMC3 Protein
Catalog #: PEH5050



Summary

Name	TMC3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7Z5M5
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NHQLKMQIQN ARSEDKKKVA QMVEARIQTQ EESTKKLPND SDLTSQLSSA HSGTPQNNGN VAHFDSGSSK SGRIETVAQS MPQSPRPGDR APSSPLPGVP KSRLEHETNR YLHGLCASTS DLHRNRSRTP MFTTTHIEDV HSEPLFRKDF QQINPPHRGP QASTLLAQGP RPHAPRYVVI NECDSYKKKH LNVWPERHFK

image



Background

Alternative Names Transmembrane channel-like protein 3

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TMC3 Protein
Catalog #: PEH5050

