


Product Name: Recombinant Human RBM23 Protein
Catalog #: PEH5060



Summary

Name	RBM23_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86U06
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KKDYPSNTTS STSNSGNETS GSSTIGETSK KKRSRSHNKS RDRKRSRSD RDRYRRRNSR SRSPGRQCRH RRSWDRRHG

image



Background

Alternative Names	Probable RNA-binding protein 23 (CAPER beta) (CAPERbeta) (RNA-binding motif protein 23) (RNA-binding region-containing protein 4) (Splicing factor SF2)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use