


Product Name: Recombinant Human ABCA12 Protein
Catalog #: PEH5062



Summary

Name	ABCAC_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86UK0
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CESETFSLIE KSCQLSDMSF GSLCEESEFD LQLLEAAELG TEIAASLLYH DNVISKKVRD LLTGDPKIN LNMDQFLEQA LQMNYLENIT QLIPIEAML HVNNSADASE KPGQLLEMFK NVEELKEDLR

image



Background

Alternative Names ATP-binding cassette sub-family A member 12 (ATP-binding cassette transporter 12) (ATP-binding cassette 12)

Note

For Research Use Only , Not for Diagnostic Use