


Product Name: Recombinant Human SHLD2 Protein
Catalog #: PEH5064



Summary

Name	SHLD2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86V20
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TEYHEIQNQC LGLFSSNAVD KSRSEAAVRK VSDLKISTDT EFLSIITSSQ VAFLAQKKDK RRSVPNKGNV NMETEPKASY GEIRIPEENS IQLDGFTEAY ESGQNQAYSL ELFSPVCPKT ENSRIHINSK KGLEEHTGSQ ELFSSSEDELP

image



Background

Alternative Names	Shieldin complex subunit 2 (Protein FAM35A) (RINN1-REV7-interacting novel NHEJ regulator 2) (Shield complex subunit 2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use