


Product Name: Recombinant Human ANKLE2 Protein
Catalog #: PEH5070



Summary

Name	ANKL2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86XL3
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSQEGLQRLE EYLTQQEIGK KAQQETGERE ASCRDKATTS GSNSISVRAF LDEDDMSLEE IKNRQNAARN NSPPTVGAFG HTRCSAFPLE QEADLIEAAE PGGPHSSRNG LCHPLNHSRT LAGKRPKAPR GEEAHLPPVS DLTVEFDKLN LQNIGRSVSK TPDESTKTKD QILTSRINAV ERDLLEPSA DQLGNNGHRT

image



Background

Alternative Names Ankyrin repeat and LEM domain-containing protein 2 (LEM domain-containing protein 4)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ANKLE2 Protein
Catalog #: PEH5070

