


Product Name: Recombinant Human TGIF2LY Protein
Catalog #: PEH5073



Summary

Name	TF2LY_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IUE0
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TQSPVEKDSP AKTQSPAQDT SIMSRNADT GRVLALPEHK KKRKGNLPAE SVKILRDWY KHRFKAYPSE EEQMLSEKT NLSLLRISNW FINARRRILP DMLQQRRNDP IGHKTGKDA HATHLQSTEA

image



Background

Alternative Names Homeobox protein TGIF2LY (TGF-beta-induced transcription factor 2-like protein) (TGFB-induced factor 2-like protein, Y-linked) (TGIF-like on the Y)

Note

For Research Use Only , Not for Diagnostic Use