


**Product Name: Recombinant Human TGIF2LX Protein**  
**Catalog #: PEH5074**



## Summary

<b>Name</b>	TF2LX_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q8IUE1
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	20 kDa
<b>Buffer</b>	Lyophilized from Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	KKRKGNLP AE SVKILRDW MY KHRFKAYP SE EEKQMLSE KT NLSLLQIS NW FINARRRIL P DMLQQRRND P IIGHKTGK DA HATHLQSTE A SVPAKSGPS G PDNVQSLPL W PLPKGQMS RE KQPDPE SAPS QKLTGIAQ PK KKVKVSVT SP SPELVSP EE

## image



## Background

**Alternative Names** Homeobox protein TGIF2LX (TGF-beta-induced transcription factor 2-like protein) (TGFB-induced factor 2-like protein, X-linked) (TGIF-like on the X)

## Note

For Research Use Only , Not for Diagnostic Use

**Product Name: Recombinant Human TGIF2LX Protein**  
**Catalog #: PEH5074**

---

