


Product Name: Recombinant Human CNGA4 Protein
Catalog #: PEH5075



Summary

Name	CNGA4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IV77
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AILQHLPERL RAEVAVSVHL STLSRVQIFQ NCEASLLEEL VLKLPQTYS PGEYVCRKGD IGQEMYIIRE GQLAVVADDG ITQYAVLGAG LYFGEISIIN IKGNMSGNRR TANIKSLGYS DLFCLSKEDL REVLSEYPQA QTIMEEKGRE ILLKMNKLDV NAEAAEIALQ EATESRLRGL DQQLDDLQTK FARLLAELES SALKIAYRIE RLEWQTREWP MPEDLAEADD EGEPEEGTSK

image



Background

Alternative Names Cyclic nucleotide-gated cation channel alpha-4 (Cyclic nucleotide-gated channel alpha-4) (CNG channel alpha-4) (CNG-4) (CNG4)

Note

Product Name: Recombinant Human CNGA4 Protein
Catalog #: PEH5075



For Research Use Only , Not for Diagnostic Use