


Product Name: Recombinant Human WDR75 Protein
Catalog #: PEH5077



Summary

Name	WDR75_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IWA0
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LNVRVMEPDP NSENIAAISQ SSVGSDLFVF KPSEPRPLYI QKGISREKVQ WGVFVPRDVP ESFTSEAYQW LNRSQFYFLT KSQSLTFST KSPEEKLTP SKQLLAEEESL PTPPFYFILG KHRQQQDEKL NETLENELVQ LPLTENIPAI SELLHTPAHV LPSAAFLCSM FVNSLLLSKE TKSKEIPED VDMEEEEKSE DSDEENDFTE KVQDTSNTGL GEDIIHQLSK

image



Background

Alternative Names	WD repeat-containing protein 75 (U3 small nucleolar RNA-associated protein 17 homolog)
--------------------------	--

Note

Product Name: Recombinant Human WDR75 Protein
Catalog #: PEH5077



For Research Use Only , Not for Diagnostic Use