


Product Name: Recombinant Human SPATA31A7 Protein
Catalog #: PEH5078



Summary

Name	S31A7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IWB4
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ELHDLMAARR SKLVHEEPRN PNCQGCKSQ RPMFPPIHKS EKSRRKPNLEK HEERLEGLRT PQLTPVRKTE DTHQDEGVQL LPSKKQPPSV SPFGENIKQI FQWIFSKKKS KPAPVTAESQ KTVKNRSRVY

image



Background

Alternative Names Spermatogenesis-associated protein 31A7 (Protein FAM75A7)

Note

For Research Use Only , Not for Diagnostic Use