

Product Name: Recombinant Human DCUN1D3 Protein
Catalog #: PEH5079



Summary

Name	DCNL3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IWE4
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GQCVTKCKNP SSTLGSKNGD REPSNKSHSR RGAGHREEQV PPCGKPGGDI LVNGTKKAEA ATEACQLPTS SGDAGRESKS NAEESLQRL EELFRRYKDE REDAILEEGM ERFCNDLCVD PTEFRVLLLA WKFQAATMCK FTRKEFFDGC KAISADSIDG ICARFPSLLT EAKQEDKFKD LYRFTFQFGL DSEEGQ RSLH REIAIALWKL VFTQNNPPVL DQWLNFLTEN PSGIGISR D TWNMFLNFTQ VIGPDLNYS EDEAWPSLFD TFVEWEMERR KREGGRGAL SSGPEGLCPE EQT

image



Background

Alternative Names	DCN1-like protein 3 (DCNL3) (DCUN1 domain-containing protein 3) (Defective in cullin neddylation protein 1-like protein 3) (Squamous cell carcinoma-related oncogene 3)
--------------------------	---

Product Name: Recombinant Human DCUN1D3 Protein
Catalog #: PEH5079



Note

For Research Use Only , Not for Diagnostic Use