


Product Name: Recombinant Human TRIM59 Protein
Catalog #: PEH5080



Summary

Name	TRI59_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IWR1
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IDDLQSAYLK EKDT PQKLE QLT DTHWTDL THLIEKLKEQ KSHSEKMIQG DKEAVLQYFK ELNDTLEQKK KSFLTALCDV GNLINQEYTP QIERMKEIRE QQLELMALTI SLQEE SPLKF

image



Background

Alternative Names	Tripartite motif-containing protein 59 (RING finger protein 104) (Tumor suppressor TSBF-1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use