


Product Name: Recombinant Human LMTK2 Protein
Catalog #: PEH5081



Summary

Name	LMTK2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IWU2
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IFDLMELNGV QADFKPATLS SSLDNPKEV ITGHFEKEKP RKIFDSEPLC LSDNLMHQDN FDPLNVQELS ENFLFLQEKV LLKGSLSSEK HINDLQTEK NAGFTEAMLE TSCRNSLDTE LQFAENKPGL SLLQENVSTK

image



Background

Alternative Names	Serine/threonine-protein kinase LMTK2 (EC 2.7.11.1) (Apoptosis-associated tyrosine kinase 2) (Brain-enriched kinase) (hBREK) (CDK5/p35-regulated kinase) (CPRK) (Kinase/phosphatase/inhibitor 2) (Lemur tyrosine kinase 2) (Serine/threonine-protein kinase KPI-2)
--------------------------	--

Note

Product Name: Recombinant Human LMTK2 Protein
Catalog #: PEH5081



For Research Use Only , Not for Diagnostic Use