


Product Name: Recombinant Human ASXL1 Protein
Catalog #: PEH5085



Summary

Name	ASXL1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IXJ9
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EIKSGLCVPG ESVRIQRGPA TRQRDGHFKK RSRPDLRTRA RRNLYKKQES EQAGVAKDAK SVASDVPLYK DGEAKTDPAG LSSPHLPGTS SAAPDLEGPE FPVESVASRI QAEPDNLARA SASPDRIPSL PQETVDQEPK DQKRKSFEQA ASASFPEKKP RLEDQRSFRN TIESVHTEKP

image



Background

Alternative Names Polycomb group protein ASXL1 (Additional sex combs-like protein 1)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ASXL1 Protein
Catalog #: PEH5085

