


Product Name: Recombinant Human SLX4 Protein
Catalog #: PEH5087



Summary

Name	SLX4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IY92
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SDVQFQTDSG EVLYAHKFVL YARCPLLIQY VNNEGFSAVE DGVLTQRVLL GDVSTEART FLHYLYTADT GLPPGLSSEL SSLAHRFGVS ELVHLCEQVP IATDSEGKPW EEKEAENCES RAENFQELLR SMWADEEEEA ETLLKSKDHE EDQENVNEAE MEEIYEFAAT QRKLLQEERA AGAGEDADWL EGGSPVSGQL

image



Background

Alternative Names	Structure-specific endonuclease subunit SLX4 (BTB/POZ domain-containing protein 12)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SLX4 Protein
Catalog #: PEH5087

