



Product Name: Recombinant Human DIS3L2 Protein
Catalog #: PEH5088



Summary

Name	DI3L2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IYB7
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SPGDKKSKNR STRGKKKSIF ETYMSKEDVS EGLKRGTLIQ GVLRLNPKKF HEAFIPSPDG DRDIFIDGVV ARNRALNGDL VVVKLLPEEH WKVVKPESND KETEAAYESD IPEELCGHHL PQQSLKSYND SPDVIVEAQF DGSDSEDGHG ITQNVLDGV KKLSVCVSEK GREDGDAPVT KDETTCISQD

image



Background

Alternative Names DIS3-like exonuclease 2 (hDIS3L2) (EC 3.1.13.-)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human DIS3L2 Protein
Catalog #: PEH5088

