


Product Name: Recombinant Human C1orf87 Protein
Catalog #: PEH5091



Summary

Name	CA087_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N0U7
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VEKPKIKEND SLKTETQTMH QKPMTDNARQ MSRDTPVPIN FTDQQTDDNP DDVKEKKHPE NNQKSENNQK LLTGANSSRF LDGNIPSQAN VHCSSVPTGD QSLSYVHGIP RRKLRDWSLE QMVRGSSDQP EDIGQSPSGT TNEDAFLLAL VRRELKSRPL SSNLLLEKLQK ELKILDPISS

image



Background

Alternative Names Uncharacterized protein C1orf87

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human C1orf87 Protein
Catalog #: PEH5091

