


Product Name: Recombinant Human NLGN4X Protein
Catalog #: PEH5092



Summary

Name	NLGNX_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N0W4
Host	E.coli
Species	Human
Molecular Mass	58 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NYGKIRGLRT PLPNEILGPV EQYLGVPYAS PPTGERRFQP PEPPSSWTGI RNTTQFAAVC PQHLDESSL HDMLPIWFTA NLDLMTYVQ DQNECLLYLN IYVPTEDDIH DQNSKKPVMV YIHGGSYMEG TGNMIDGSIL ASYGNVIVIT INYRLGILGF LSTGDQAAKG NYGLLDQIQA LRWIEENVGA FGGDPKRVTI FGSGAGASCV SLLTSLHYSE GLFQKAIQS GTALSSWAVN YQPAKYTRIL ADKVGCNMLD TDMVECLRN KNYKELIQQT ITPATYHIAF GPVIDGDVIP DDPQILMEQG EFLNYDIMLG VNQGEGLKFV DGIVDNEDGV TPNDFDFSVS NFVDNLYGYP EGKDTLRETI KFMYPDWADK ENPETRRKTL VALFTDHQWV APAVATADLH AQYGSPTYFY AFYHHCQSEM KPSWADSAHG DEVPYVFGIP MIGPTLFC NFSKNDVMLS AVVMYWTNF AKTGDPNQPV PQDTKFIHTK

image



Background

Product Name: Recombinant Human NLGN4X Protein
Catalog #: PEH5092



Alternative Names

Neurologin-4, X-linked (Neurologin X) (HNLX)

Note

For Research Use Only , Not for Diagnostic Use