


Product Name: Recombinant Human VPS52 Protein
Catalog #: PEH5096



Summary

Name	VPS52_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N1B4
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ARELVL RAGT SDMEEEEGPL AGGPGLQEPL QLGELDITSD EFILDEV DVH IQANLEDELV KEALKTGV DL RHYSKQVELE LQQIEQKSIR DYIQESENIA SLHNQITACD AVLERMEQML GAFQSDLSSI SSEIRTLQEQ SGAMNIRLRN RQAVRGKLG E LVDGLVPSA LVTAILEAPV TEPRFLEQLQ ELDAKAAAVR

image



Background

Alternative Names Vacuolar protein sorting-associated protein 52 homolog (SAC2 suppressor of actin mutations 2-like protein)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human VPS52 Protein
Catalog #: PEH5096

