


Product Name: Recombinant Human REXO1 Protein
Catalog #: PEH5098



Summary

Name	REXO1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N1G1
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ALEKEVDFDS DPMEECLRIF NESTSVKTED RGRLARQPPK EEKSEEKGLS GLTTLFPGQK RRISHLSKQG QEVEPPRRGP

image



Background

Alternative Names	RNA exonuclease 1 homolog (EC 3.1.-.-) (Elongin-A-binding protein 1) (EloA-BP1) (Transcription elongation factor B polypeptide 3-binding protein 1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use