


Product Name: Recombinant Human ZNF710 Protein
Catalog #: PEH5100



Summary

Name	ZN710_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N1W2
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PPVRLVPKVK FEKVEEEQE VYEVSVPGDD KDAGPAEAPA EAASGGCDAL VQSSAVKMID LSAFSRKPRP LRHLPRTPRP ELNVAPYDPH FPAPARDGFP EPSMALPGPE ALPTECGFEP PHLAPLSDPE APSMESPEPV KPEQGFVWQE ASEFEADTAG

image



Background

Alternative Names Zinc finger protein 710

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZNF710 Protein
Catalog #: PEH5100

