

Summary

Name	CLSR_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N2M8
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KLRGMMVDYK KRAERRREYY EKIKKDP AQF LQVHGRACKV HLDSAVALAA ESPVNMMMPWQ GDTNNMIDRF DVRAHLDHIP DYTPPLTTI SPEQESDERK CNYERYRGLV QNDFAGISEE QCLYQIYIDE LYGGLQRPSE DEKKKLAEEK ASIGYTYEDS TVAKVEKAAE KPEEEESAAE EESNSDEDEV IPDIDVEVDV DELNQEQVAD LNKQATTYGM ADGDFVRMLR KDKEEA EAIK HAKALEEEKA MYSGRRSRRQ

image



Background

Alternative Names	CLK4-associating serine/arginine rich protein (Splicing factor, arginine/serine-rich 16) (Suppressor of white-apricot homolog 2)
--------------------------	--

Product Name: Recombinant Human CLSRP Protein
Catalog #: PEH5102



Note

For Research Use Only , Not for Diagnostic Use