


Product Name: Recombinant Human ZNF276 Protein
Catalog #: PEH5110



Summary

Name	ZN276_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N554
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LCRAPGQLGE KQLPSSTSD RVKDEFSDLS EGDVLESDEN DKKQNAQSSD ESFEPYPERK VSGKKSEKSE AKKSEEPRIK KKP GPKPGWK KKLRCEREEL PTIYKCPYQG CTAVYRGADG MKKHIKEHHE EVRERPCPHP

image



Background

Alternative Names Zinc finger protein 276 (Zfp-276) (Zinc finger protein 477)

Note

For Research Use Only , Not for Diagnostic Use