


Product Name: Recombinant Human ZGPAT Protein
Catalog #: PEH5111



Summary

Name	ZGPAT_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N5A5
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GDLKELIELT EASLVSVRKS SLLAALDEER PGRQEDAEYQ AFREAITEAV EAPAAARGSG SETVPKAEAG

image



Background

Alternative Names	Zinc finger CCCH-type with G patch domain-containing protein (G patch domain-containing protein 6) (Zinc finger CCCH domain-containing protein 9) (Zinc finger and G patch domain-containing protein)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use