


Product Name: Recombinant Human MSL3 Protein
Catalog #: PEH5115



Summary

Name	MS3L1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N5Y2
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FHSGEKVLCF EPDPTKARVL YDAKIVDVIV GKDEKGRKIP EYLIHFNGWN RSWDRWAAED HVLRTDENR RLQRKLARKA VARLRSTGRK KKRCRLPGVD SVLKGLPTEE KDENDENSLS SSSDCSENKD EEISEESDIE

image



Background

Alternative Names	Male-specific lethal 3 homolog (Male-specific lethal-3 homolog 1) (Male-specific lethal-3 protein-like 1) (MSL3-like 1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use