


Product Name: Recombinant Human ATP11C Protein
Catalog #: PEH5127



Summary

Name	AT11C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NB49
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EIEEYELLHT LNFDAVRRRM SVIVKTQEGD ILLFCKGADS AVFPRVQNHE IELTKVHVER NAMDGYRTL C VAFKEIAPDD YERINRQLIE AKMALQDREE KMEKVFDDIE TNMNLIGATA VEDKLQDQAA ETIEALHAAG LKVWVLTGDK METAKSTCYA CRLFQTNTL ELTTKTIEE

image



Background

Alternative Names	Phospholipid-transporting ATPase IG (EC 7.6.2.1) (ATPase IQ) (ATPase class VI type 11C) (P4-ATPase flippase complex alpha subunit ATP11C)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ATP11C Protein
Catalog #: PEH5127

