


Product Name: Recombinant Human GOLM1 Protein
Catalog #: PEH5128



Summary

Name	GOLM1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NBJ4
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DLQTRIMELE GRVRRAAAER GAVELKKNEF QGELEKQREQ LDKIQSSHNF QLESVNKLYQ DEKAVLVNNI TTGERLIRVL QDQLKTLQRN YGRLQQDVLQ FQKNQTNLER KFSYDLSQCI NQMKEVKEQC EERIEEVTKK GNEAVASRDL SENNDQRQQL QALSEPQPRL QAAGLPHEV PQGKGNVLGN SKSQTPAPSS EVLVLSKRQV EKEETNEIQV VNEEPQRDRL PQEPGREQVV

image



Background

Alternative Names Golgi membrane protein 1 (Golgi membrane protein GP73) (Golgi phosphoprotein 2)

Note

Product Name: Recombinant Human GOLM1 Protein
Catalog #: PEH5128



For Research Use Only , Not for Diagnostic Use