


Product Name: Recombinant Human HJURP Protein
Catalog #: PEH5129



Summary

Name	HJURP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NCD3
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GHGRNRYDEI KEEFDKLHQK YCLKSPGQMT VPLCIGVSTD KASMEVRYQT EGFLGKLNPDPHFQGFQKLPSSPLGCRKSL LGSTAIEAPS STCVARAITR DGTRDHQFPA KRPRLSEPQG SGRQGNSLGA SDGVDNTVRP GDQGSSSQPN SEERGENTSY

image



Background

Alternative Names	Holliday junction recognition protein (14-3-3-associated AKT substrate) (Fetal liver-expressing gene 1 protein) (Up-regulated in lung cancer 9)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human HJURP Protein
Catalog #: PEH5129

