


Product Name: Recombinant Human DYNC2H1 Protein
Catalog #: PEH5133



Summary

Name	DYHC2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NCM8
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LKRQSHLQAG VSKLNEAKAL VDELNRKAGE QSVLLKTKQD EADAALQMIT VSMQDASEQK TELERLKHRI AEEVVKIEER KNKIDDELKE VQPLVNEAKL AVGNIKPESL SEIRSLRMPP DVIRDILEGV

image



Background

Alternative Names	Cytoplasmic dynein 2 heavy chain 1 (Cytoplasmic dynein 2 heavy chain) (Dynein cytoplasmic heavy chain 2) (Dynein heavy chain 11) (hDHC11) (Dynein heavy chain isotype 1B)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human DYNC2H1 Protein
Catalog #: PEH5133

