


Product Name: Recombinant Human CNNM3 Protein
Catalog #: PEH5135



Summary

Name	CNNM3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NE01
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LED CFMLDAS TVLDFGVLAS IMQSGHTRIP VYEEERSNIV DMLYLKDLAF VDPEDCTPLS TITRFYNHPL HFVFNDTKLD AVLEEFKRGK SHLAIVQKVN NEGEGDPFYE VLGLVTLEDV IEEIIRSEIL DESEDYRDTV VKRKPASLMA PLKRKEEFS LFKVSDDEYKV TISPQLLLAT

image



Background

Alternative Names	Metal transporter CNNM3 (Ancient conserved domain-containing protein 3) (Cyclin-M3)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CNNM3 Protein
Catalog #: PEH5135

