


Product Name: Recombinant Human RHOXF1 Protein
Catalog #: PEH5143



Summary

Name	RHXF1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NHV9
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ISPTPQLGAA SSAEGHVGQG APGLMGNMNP EGGVNHENGM NRDGGMPIPEG GGGNQEPRQQ PQPPPEEPAQ AAMEGPQPEN MQPRTRRTKF TLLQVEELES VFRHTQYPDV PTRRELAENL

image



Background

Alternative Names Rhox homeobox family member 1 (Ovary-, testis- and epididymis-expressed gene protein) (Paired-like homeobox protein PEPP-1)

Note

For Research Use Only , Not for Diagnostic Use