


Product Name: Recombinant Human PCDH19 Protein
Catalog #: PEH5146



Summary

Name	PCD19_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TAB3
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GNSLKDSGHE ESDQTDSEHD VQRSLYCDTA VNDVLNTSVT SMGSQMPDHD QNEGFHCREE CRILGHSDRC WMPRNPMPPIR SKSPEHVRNI IALSIEATAA DVEAYDDCGP TKRTFATFGK DVSDHPAEER PTLKGKRTVD VTICSPKVNS VIREAGNGCE AISPVT SPLH LKSSLPTKPS VSYTIALAPP ARDLEQYVNN VNNGPTRPSE AEPRGADSEK VMHEVSPILK EGRNKESPGV

image



Background

Alternative Names Protocadherin-19

Note

Product Name: Recombinant Human PCDH19 Protein
Catalog #: PEH5146



For Research Use Only , Not for Diagnostic Use