


Product Name: Recombinant Human GABPB2 Protein
Catalog #: PEH5148



Summary

Name	GABP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TAK5
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ANPVTDPVSM AAPFIFTSGE VVNLASLISS TNTKTTSGDP HASTVQFSNS TTSVLATLAA LAEASVPLSN SHRATANTEE IIEGNSVDSS IQQVMGSGGQ RVITIVTDGV PLGNIQTSIP TGGIGQPFIV TVQDGQQVLT VPAGKVAEET VIKEEEEEKL PLTKKPRIGE KTNSVEESKE GNERELLQQQ LQEANRRAQE YRHQLLKEQ EAEQYRLKLE AIARQQPNGV DFTMVEEVAE

image



Background

Alternative Names GA-binding protein subunit beta-2 (GABP subunit beta-2) (GABPB-2)

Note

Product Name: Recombinant Human GABPB2 Protein
Catalog #: PEH5148



For Research Use Only , Not for Diagnostic Use