


Product Name: Recombinant Human GOLGA5 Protein
Catalog #: PEH5151



Summary

Name	GOGA5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TBA6
Host	E.coli
Species	Human
Molecular Mass	64 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ELLFDFLNSS QKEPTGRVEI RKEKGKTPVF QSSQTSSVSS VNPSVTTIKT IEENSFGSQT HEAASNSDSS HEGQEESKE NVSSNAACPD HTPTPNDDGK SHELNLRLLE NQLLRNEVQS LNQEMASLLQ RSKETQEELN KARARVEKWN ADHSKSDRMT RGLRAQVDDL TEAVAAKDSQ LAVLKVRLQE ADQLLSTRTE ALEALQSEKS RIMQDQSEGN SLQNQALQTF QERLHEADAT LKREQESYKQ MQSEFAARLN KVEMERQNLA EAITLAERKY SDEKKRVDEL QQQVKLYKLN LESSKQELID YKQKATRILQ SKEKLINSLK EGSGFEG LDS STASSMELEE LRHEKEMQRE EIQLKMGQIH QLRSELQDME AQQVNEAESA REQLQDLHDQ IAGQKASKQE LETELERLKQ EFHYIEEDLY RTKNTLQSRI KDRDEEIQKL RNQLTNKTLS NSSQSELENR LHQLTETLIQ KQTMLESLSST EKNSLVFQLE RLEQQMNSAS GSSSNGSSIN MSGIDNGEGT RLRNVPVLFN

image



Background

Product Name: Recombinant Human GOLGA5 Protein
Catalog #: PEH5151



Alternative Names

Golgin subfamily A member 5 (Cell proliferation-inducing gene 31 protein)
(Golgin-84) (Protein Ret-II) (RET-fused gene 5 protein)

Note

For Research Use Only , Not for Diagnostic Use