


Product Name: Recombinant Human ZCRB1 Protein
Catalog #: PEH5153



Summary

Name	ZCRB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TBF4
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RAAEFIRRRN YFDKSKCYEC GESGHLASYAC PKNMLGEREP PKKKEKKKKK KAPEPEEEIE EVEESEDEGE DPALDSLSQA IAFQQAKIEE EQKKWKPSG VPSTSDDSRR PRIKKSTYFS

image



Background

Alternative Names	Zinc finger CCHC-type and RNA-binding motif-containing protein 1 (U11/U12 small nuclear ribonucleoprotein 31 kDa protein) (U11/U12 snRNP 31 kDa protein) (U11/U12-31K)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZCRB1 Protein
Catalog #: PEH5153

