


Product Name: Recombinant Human CCDC110 Protein
Catalog #: PEH5154



Summary

Name	CC110_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TBZ0
Host	E.coli
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DEVSVLLSA SKILNSSEGV KESGCSDEY GCIAEENQI QPQSALKVLQ QQLESFQALR MQTLQNVSMV QSEISEILNK SIIEVENPQF SSEKNLVFGT RIEKDLPTEN QEENLSMEKS HFFEDSKTLH SVEEKLSGDS VNSLPQSVNV PSQIHSEDTL TLRTSTDNLS SNIIHPSEN SDILKNYNNF YRFLPTAPPN VMSQADTVIL DSKITVPFL KHGFCENLDD ICHSIKQMKE ELQKSHDGEV ALTNELQTLQ TDPDVHRNGK YDMSPIHQDK MNFIKEENLD GNLNEDIKSK RISELEALVK KLLPFRETVS KFHVHFCRKC KKLSKSEMHR GKKNEKNNKE IPITGKNITD

image



Background

Alternative Names Coiled-coil domain-containing protein 110 (Cancer/testis antigen 52) (CT52)

Product Name: Recombinant Human CCDC110 Protein
Catalog #: PEH5154



(Cancer/testis antigen KM-HN-1)

Note

For Research Use Only , Not for Diagnostic Use