


Product Name: Recombinant Human MGARP Protein
Catalog #: PEH5160



Summary

Name	HUMMR_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TDB4
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QAKHTEHKTN LKEKTKAEIH PFQGEKENVA ETEKASSEAP EELIVEAEVV DAEESPSATV VVIKEASACP GHVEAAPETT AVSAETGPEV TDAAARETTE VNPETTPEVT NAALDEAVTI DNDKDTTKNE TSDEYAELEE ENSPAECESS AGDDLQEEAS VGSEAASAQG

image



Background

Alternative Names	Protein MGARP (Corneal endothelium-specific protein 1) (CESP-1) (Hypoxia up-regulated mitochondrial movement regulator protein) (Mitochondria-localized glutamic acid-rich protein) (Ovary-specific acidic protein)
--------------------------	---

Note

Product Name: Recombinant Human MGARP Protein
Catalog #: PEH5160



For Research Use Only , Not for Diagnostic Use