


Product Name: Recombinant Human FAT3 Protein
Catalog #: PEH5161



Summary

Name	FAT3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TDW7
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KIVDGLGGEH QEMTTFHPES PRILTARRGV VVCSVAPNLP AVSPCRSDCD SIRKNGWDAG TENKGVDDPG EVTCFAGSNK GSNSEVQSLS SFQSDSGDDN AYHWDTSWDM PGARLSDIEE VPNYENQDGG SAHQGSTREL ESDYYLGGYD IDSEYPPPH

image



Background

Alternative Names	Protocadherin Fat 3 (hFat3) (Cadherin family member 15) (FAT tumor suppressor homolog 3)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human FAT3 Protein
Catalog #: PEH5161

