


Product Name: Recombinant Human CSGALNACT1 Protein
Catalog #: PEH5162



Summary

Name	CGAT1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TDX6
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EQLALPRANS PTGKEGYQAV LQEWEEQHRN YVSSLKRQIA QLKEELQERS EQLRNGQYQA SDAAGLGLDR SPPEKTQADL LAFLHSQVDK AEVNAGVKLA TEYAAVPFDS FTLQKVYQLE TGLTRHPEEK PVRKDKRDEL VEAIESALET LNSPAENSPN HRPYTASDFI

image



Background

Alternative Names	Chondroitin sulfate N-acetylgalactosaminyltransferase 1 (CsGalNAcT-1) (EC 2.4.1.174) (Chondroitin beta-1,4-N-acetylgalactosaminyltransferase 1) (Beta4GalNAcT-1)
--------------------------	--

Note



Product Name: Recombinant Human CSGALNACT1 Protein
Catalog #: PEH5162

For Research Use Only , Not for Diagnostic Use