


Product Name: Recombinant Human LMX1A Protein
Catalog #: PEH5163



Summary

Name	LMX1A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TE12
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EESLCKSAHG AGKGTAEEGK DHKRPKRPT ILTTQQRRAF KASFEVSSKP CRKVRETLAA ETGLSVRVVQ VWFQNQRAKM KKLARRQQQQ

image



Background

Alternative Names LIM homeobox transcription factor 1-alpha (LIM/homeobox protein 1.1) (LMX-1.1) (LIM/homeobox protein LMX1A)

Note

For Research Use Only , Not for Diagnostic Use