


Product Name: Recombinant Human FBF1 Protein
Catalog #: PEH5165



Summary

Name	FBF1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TES7
Host	E.coli
Species	Human
Molecular Mass	63 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LSATEPATCF PSTQKPTEPS VPVQPLLPE LARSLLPSTE YQKQLLAAQV QLQCSPAELQ AELLHSQARL AELEAQVRKL ELERAQHELL LGSLLQQHQA DLELIESAHR SRIKVLETSY QQREERLRRE NEELSARYLS QCQEAQARA ELTAQHQRRL AAIAQEKDQE MERLRELQRA SILDMMRDHE EQLQRLKLLK DREVDAA TSA TSHTRSLNSI IHQMEKFSS LHELSSRVEA SHLTTSQERE LGIRQRDEQL RALQERLGQQ QRDMEEERSR QQEVIGKMEA RLNEQSRLLE QERWRVTAEQ SKAESMQRAL EEQRKVTAQQ MAMERAELER AKSALLEEQK SVMLKCGEER RRLAAEWA EF SAQQKLSKER AEREERALQ VDTQREGTLI SLAKQAE LKI RASELR AEEK QLAAERAALE QERQELRLEK ERINATALRV KLRAEEVESM SKVASEKYEE GERALREAQQ VQAEQQARLQ AVQQQQERLR KQEQHMHQEH LSLAQQLQL

image



Background

Product Name: Recombinant Human FBF1 Protein
Catalog #: PEH5165



Alternative Names

Fas-binding factor 1 (FBF-1) (Protein albatross)

Note

For Research Use Only , Not for Diagnostic Use