


Product Name: Recombinant Human PNISR Protein
Catalog #: PEH5166



Summary

Name	PNISR_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TF01
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SDPEMTEEEK EYQMMLLTKM LLTEILLDVT DEEIYYVAKD AHRKATKAPA KQLAQSSALA SLTGLGGLGG YGSGDSEDER SDRGSESSDT DDEELRHRIR QKQEAFWRKE KEQQLLHDKQ MEEKQQTER VTKEMNEFIH KEQNSLSLLE AREADGDVVN EKKRTPNETT SVLEPKKEHK

image



Background

Alternative Names	Arginine/serine-rich protein PNISR (PNN-interacting serine/arginine-rich protein) (SR-related protein) (SR-rich protein) (Serine/arginine-rich-splicing regulatory protein 130) (SRrp130) (Splicing factor, arginine/serine-rich 130) (Splicing factor, arginine/serine-rich 18)
--------------------------	--

Product Name: Recombinant Human PNISR Protein
Catalog #: PEH5166



Note

For Research Use Only , Not for Diagnostic Use