


Product Name: Recombinant Human ATP8B4 Protein
Catalog #: PEH5168



Summary

Name	AT8B4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TF62
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	WHKMLEDANA ATEERDERIA GLYEEIERDL MLLGATAVED KLQEGVIETV TSLSLANIKI WVLGTGDKQET AINIGYACNM LTDDMNDVFV IAGNNAVEVR EELRKAKQNL FGQNRNFSNG HVVCEKKQQL ELDSIVEETI TGDYALIING HSLAHALESD

image



Background

Alternative Names	Probable phospholipid-transporting ATPase IM (EC 7.6.2.1) (ATPase class I type 8B member 4) (P4-ATPase flippase complex alpha subunit ATP8B4)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ATP8B4 Protein
Catalog #: PEH5168

