


Product Name: Recombinant Human SLC20A1 Protein
Catalog #: PEH5171



Summary

Name	S20A1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WUM9
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EIKCSPESP LMEKKNLKE DHEETKLSVG DIENKHPVSE VGPATVPLQA VVEERTVSFK LGDLEEAPER ERLPSVDLKE ETSIDSTVNG AVQLPNGNLV QFSQAVSNQI NSSGHYQYHT VHKDSGLYKE

image



Background

Alternative Names	Sodium-dependent phosphate transporter 1 (Gibbon ape leukemia virus receptor 1) (GLVR-1) (Leukemia virus receptor 1 homolog) (Phosphate transporter 1) (PiT-1) (Solute carrier family 20 member 1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SLC20A1 Protein
Catalog #: PEH5171

