

Product Name: Recombinant Human LRRN4 Protein
Catalog #: PEH5172



Summary

Name	LRRN4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WUT4
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LSRAADTMCA PAAGSSGPFS ASLSLSQLPG VCQSDQSTTL GASHPPCFNR STYAQGTTVA PSAAPATRPA GDQQSVSKAP NVGSRTIAAW PHSDAREGTA PSTTNSVAGH SNSSVFPRAA STTRTQHRGE HAPELVLEPD

image



Background

Alternative Names Leucine-rich repeat neuronal protein 4 (Neuronal leucine-rich repeat protein 4) (NLRR-4)

Note

For Research Use Only , Not for Diagnostic Use