


Product Name: Recombinant Human STAG1 Protein
Catalog #: PEH5175



Summary

Name	STAG1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WVM7
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LAKRRGRNSP NGNLIRMLVL FFLESELHEH AAYLVDSLWE SSQELLKDWE CMTELLLEEP VQGEEAMSDR QESALIELMV CTIRQAEEAH PPVGRGTGKR VLTAKERKTQ IDDRNKLTEH FIITLPMLLS KYSADA EKVA NLLQIPQYFD LEIYSTGRME KHL DALLKQI K FVVEKHVES DVLEACSKTY SILCSEEYTI QNRVDIARSQ LIDEFVDRFN HSVEDLLQEG

image



Background

Alternative Names	Cohesin subunit SA-1 (SCC3 homolog 1) (Stromal antigen 1)
--------------------------	---

Note

Product Name: Recombinant Human STAG1 Protein
Catalog #: PEH5175



For Research Use Only , Not for Diagnostic Use