


Product Name: Recombinant Human ZFPM2 Protein
Catalog #: PEH5176



Summary

| | |
|------------------------------|--|
| Name | FOG2_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q8WW38 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 10 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | HQELHVPSGK LPRESMEHS PSATEDSLQP ATDLLTRSEL PQSQKAMQTK DASSDTELDK CEKKTQLFLT |

image



Background

| | |
|--------------------------|--|
| Alternative Names | Zinc finger protein ZFPM2 (Friend of GATA protein 2) (FOG-2) (Friend of GATA 2) (hFOG-2) (Zinc finger protein 89B) (Zinc finger protein multitype 2) |
|--------------------------|--|

Note

For Research Use Only , Not for Diagnostic Use