


**Product Name: Recombinant Human LMO7 Protein**  
**Catalog #: PEH5178**



## Summary

<b>Name</b>	LMO7_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q8WWI1
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	52 kDa
<b>Buffer</b>	Lyophilized from Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	VSRSDQFSD MRISINQTPG KSLDFGFTIK WDIPGIFVAS VEAGSPAEF QLQVDDEIIA INNTKFSYND SKEWEEAMAK AQETGHLVMD VRRYGKAGSP ETKWIDATSG IYNSEKSSNL SVTTDFSESL QSSNIESKEI NGIHDESNF ESKASESISL KNLKRRSQFF EQGSSDSVVP DLPVPTISAP SRWVWDQEEE RKRQERWQKE QDRLLQEKYQ REQEKLREEW QRAKQEAERE NSKYLDEELM VLSSNSMSLT TREPSTATWE ATWSEGSKSS DREGTRAGEE ERRQPQEEVV HEDQGKKPQD QLVIERERKW EQQLQEEQEQ KRLQAEAEQ KRPAAEQKRQ AEIERETSVR IYQYRRPVDS YDIPKTEEAS SGFLPGDRNK SRSTTELDY STNKNGNKKY LDQIGNMTSS QRRSKEQVP SGAELERQQI

## image



## Background

**Product Name: Recombinant Human LMO7 Protein**  
**Catalog #: PEH5178**



---

**Alternative Names**

LIM domain only protein 7 (LMO-7) (F-box only protein 20) (LOMP)

**Note**

For Research Use Only , Not for Diagnostic Use