


Product Name: Recombinant Human PRPF31 Protein
Catalog #: PEH5179



Summary

Name	PRP31_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WWY3
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GSYGEEEEEP AIEDVQEETQ LDLSGDSVKT IAKLWDSKMF AEIMMKIEEY ISKQAKASEV MGPVEAAPEY RVIVDANNLT VEIENELNII HKFIRDKYSK RFPELESLVP NALDYIRTVK ELGNSLDKCK NNENLQQILT

image



Background

Alternative Names	U4/U6 small nuclear ribonucleoprotein Prp31 (Pre-mRNA-processing factor 31) (Serologically defined breast cancer antigen NY-BR-99) (U4/U6 snRNP 61 kDa protein) (Protein 61K) (hPrp31)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PRPF31 Protein
Catalog #: PEH5179

