


Summary

Name	RUFY2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WXA3
Host	E.coli
Species	Human
Molecular Mass	50 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LDSQVGVIDF SMYLNKNEEDI GNKERNVQIA AILDQKNYVE ELNRQLNSTV SSLHSRVDSL EKSNTKLIIE LAIAKNNIIK LQEENHQLRS ENKLILMKTQ QHLEVTKVDV ETELQTYKHS RQGLDEMYNE ARRQLRDESQ LRQDVENELA VQVSMKHEIE LAMKLEKDI HEKQDTLIGL RQQLLEEVKAI NIEMYQKLQG SEDGLKEKNE IARLEEKTN KITAAMRQLE QRLQQAQEAQ MEAEDEDEKY LQECLSKSDS LQKQISQKEK QLVQLETDLK IEKEWRQTLQ EDLQKEKDAL SHLRNETQQI ISLKKEFLNL QDENQQLKKI YHEQEQALQE LGNKLSESKL KIEDIKEANK ALQGLVWLKD KEATHCKLCE KEFSLSKRKH HCRNCGEIFC NACSDNELPL PSSPKPVRVC

image



Background

Product Name: Recombinant Human RUFY2 Protein
Catalog #: PEH5181



Alternative Names

RUN and FYVE domain-containing protein 2 (Rab4-interacting protein related)

Note

For Research Use Only , Not for Diagnostic Use